

## Infinitely better compatibility with any LC

ZORBAX RRHD columns are available in four C18 phases for the ultimate flexibility in fine-tuning your selectivity to meet the needs of your analysis. Eclipse Plus C18 is a good place to start for most applications, offering high performance and excellent peak shapes with acids, bases and neutrals over the pH range 2 - 9. Try Stablebond C18 (non-end-capped) for alternate selectivity and low-pH (between pH 1 - 2) stability. The Eclipse XDB C18 phase features an extra dense bonding for long lifetime and performance over a wide pH range. Try Extend-C18 for high pH applications. You will achieve the same selectivity using ZORBAX RRHD columns as you would on ZORBAX Rapid Resolution High Throughput (RRHT, 1.8  $\mu$ m), 3.5, and 5  $\mu$ m columns with the same bonded phase so you are assured of method scalability.

### ZORBAX Rapid Resolution High Definitions (RRHD) Columns for High Pressure Use (Maximum Pressure: 1200 bar)

	Eclipse Plus C18	Eclipse Plus C8	StableBond-C18	StableBond-C8	Eclipse XDB-C18	Extend-C18
ZORBAX RRHD 2.1 x 150 mm, 1.8 $\mu$ m	959759-902	959759-906	859700-902	859700-906	981759-902	759700-902
ZORBAX RRHD 2.1 x 100 mm, 1.8 $\mu$ m	959758-902	959758-906	858700-902	858700-906	981758-902	758700-902
ZORBAX RRHD 2.1 x 50 mm, 1.8 $\mu$ m	959757-902	959757-906	857700-902	857700-906	981757-902	757700-902
ZORBAX RRHD 3.0 x 150 mm, 1.8 $\mu$ m	959759-302	959759-306	859700-302	859700-306	981759-302	759700-302
ZORBAX RRHD 3.0 x 100 mm, 1.8 $\mu$ m	959758-302	959758-306	858700-302	858700-306	981758-302	758700-302
ZORBAX RRHD 3.0 x 50 mm, 1.8 $\mu$ m	959757-302	959757-306	857700-302	857700-306	981757-302	757700-302



**Agilent Technologies**