

Millions of Chromacol vials have been used with this instrument for many years in its characteristic black cassette system. Our MICRO+™ vials with a capacities of 300µL and 900µL are both easy to use and very safe for the autosampler needle. Our Instrument Select product code HPL includes Chromacol 2-CV vials and 11-AC7 caps.

Key to products ■ Vial ■ Cap ■ Combination Pack

Bold items in grey are our standard products

PART NUMBER	DESCRIPTION	12 x 32mm crimp and snap cap vials and caps
HPL	2mL crimp top vial with a type 7 rubber/PTFE seal	
2-CV	2mL, clear glass, crimp top, flat bottom, wide neck vial. With 'write-on' patch.	
1.1-CTVG	1.1mL, clear glass, crimp top, tapered, wide neck vial. Made from inert CHROMACOL GOLD™ grade glass. Needs PTFE support, reference TTS-313. See page 48.	
09-FIV	900µL Micro+™ crimp top, fused insert vial, clear.	
06-PECV	600µL, polyethylene vial with an internal tapered profile.	
03-FIV	300µL Micro+™, one piece, clear glass insert and vial combination. With 'write-on' patch.	
11-AC7	Aluminium crimp cap fitted with a type 7, rubber/PTFE seal. Shore hardness value 60. Other colours and seals available see page 34/35 for details.	
2-RV	2mL snap cap vial with write on patch - clear	
11-PSN(B)-8RTI	11mm polyethylene snap cap pre-fitted with a red rubber/PTFE seal. Shore hardness value 58.	

PART NUMBER	DESCRIPTION	6 x 32mm crimp top vial and cap
03-CVG	300µL, clear glass, crimp top, round bottom vial, made from inert CHROMACOL GOLD™ grade glass needs SV-S4 support. See page 30/31.	
8-AC7	Aluminium crimp cap fitted with a type 7, rubber/PTFE seal. Shore hardness value 60. Other colours and seals available see page 32/33 for details.	

PART NUMBER	DESCRIPTION	12 x 32mm wide neck screw top vials and cap
2-SVW	2mL, clear glass, flat bottom, screw top, wide neck vial. With 'write-on' patch.	
09-FISV	900µL Micro+™ screw top, fused insert vial, clear.	
03-FISV	300µL Micro+™, one piece, clear glass insert and vial combination. With 'write-on' patch.	
9-SC(B)-8RTI	9mm screw cap (blue) prefitted rubber PTFE seal. Shore hardness value 58. Other colours and seals available see page 38/39 for details.	
Note 4-SV	4mL wash vial, 15 x 46mm, clear glass, screw top, flat bottom. See page 40.	

