

CTC LC PAL/GC PAL

As a result of extensive testing, the recommended seals for these units are our ST-I44 and ST-I01 formulations. These are blue silicone elastomers with a PTFE coating. CTC instruments are compatible with our WebSeal™ 96 well titer plate products. Instrument Select is only available for the I05 tray version and uses product code CTCG and CTCL. A detailed brochure on CTC Consumables is also available.

Key to products ■ Vial ■ Cap ■ Seal ■ Combination Pack

Bold items in grey are our standard products



| PART NUMBER | DESCRIPTION | 7 x 40mm crimp top vials and caps |
|----------------------|--|-----------------------------------|
| 08-CPV(A) | 800µL, amber glass, crimp top, flat bottomed vial. | |
| 07-CPV | 700µL, clear glass, crimp top, tapered vial. | |
| 8-AC(B)-ST144 | Aluminium crimp cap (blue) fitted with a blue silicone/red PTFE seal. Other colours and seals available see page 32/33 for details. | |

| PART NUMBER | DESCRIPTION | 22 x 45mm crimp top vial cap and seal |
|---------------|--|---------------------------------------|
| 10-CV | 10mL, clear glass, crimp top vial. | |
| 20-MCB | Tin plate magnetic crimp cap, no seal. | |
| 20-ST3 | Blue silicone/PTFE, 3mm thick seal. Shore hardness value 45. | |

| PART NUMBER | DESCRIPTION | 12 x 32mm crimp top vials and cap |
|-----------------------|---|-----------------------------------|
| CTCG | 2mL crimp top vial with a blue silicone/red PTFE seal | |
| CTCL | 2mL crimp top vial with a blue silicone/red PTFE seal | |
| 2-CV | 2mL, clear glass, crimp top, flat bottom, wide neck vial. With 'write-on' patch. | |
| 1.1-CTVG | 1.1mL, clear glass, crimp top, tapered, wide neck vial. Made from inert CHROMACOL GOLD™ grade glass. Needs TTS-312 PTFE support. See page 48. | |
| 09-FIV | 900µL Micro+™ crimp top, fused insert vial, clear. | |
| 03-FIV | 300µL Micro+™, one piece, clear glass insert and vial combination. With 'write-on' patch. | |
| 11-AC(B)-ST144 | Aluminium crimp cap, blue, fitted with a blue silicone/red PTFE seal. | |

| PART NUMBER | DESCRIPTION | 12 x 32mm wide neck screw top vials and cap |
|----------------------|---|---|
| 2-SVW | 2mL, clear glass, flat bottom, screw top, wide neck vial. With 'write-on' patch. | |
| 09-FISV | 900µL Micro+™ screw top, fixed insert vial, clear. | |
| 03-FISV | 300µL Micro+™, one piece, clear glass insert and vial combination. With 'write-on' patch. | |
| 9-SC(B)-ST101 | Screw cap, blue, pre-fitted with a blue silicone/PTFE seal. Shore hardness value 30. Other colours and seals available see page 38/39 for details. | |

| PART NUMBER | DESCRIPTION | 6 x 32mm crimp top vial and cap |
|----------------------|--|---------------------------------|
| 03-CVG | 300µL, clear glass, crimp top, round bottom vial, made from inert CHROMACOL GOLD™ grade glass. Uses the SV-TSP support sleeve. See page 30/31. | |
| 8-AC(B)-ST144 | Aluminium crimp cap, blue, fitted with a blue silicone/red PTFE seal. | |

