

Tips with a Perfect Fit.

4060-1002 Example catalog number

Use this tip chart to determine which tips are most appropriate for use with your Ovation and MLA pipettes.

Ovation ESC (electronic single channel) and Ovation E8 & E12 (multichannel):		0.5-20µL	2-125µL 5-250µL	25-850 25-1250µL		
Ovation Macro5 & Macro10 (macro volume):					100-5000µL	0.2-10mL
Ovation QS (quick-set adjustable volume) and Ovation M (mechanical volume set):		0.2-2µL 1.0-10µL 2.0-20µL	10-100µL 20-200µL	100-1000µL		
Ovation F1 & Ovation F2 (fixed volume):		1µL to 20µL	21µL to 200µL	201µL to 1000µL		
All MLA models:			2µL to 200µL	201µL to 1000µL	1mL to 5mL	6mL to 10mL
TIP SIZE:		MICRO	SMALL	LARGE	5mL³	10mL
NON-STERILE	VistaRak 192 tips/rack, 5 racks	4060-1002	4060-2004	4060-3004		
	VistaRak, Low Retention 192 tips/rack, 5 racks	4070-1002LR	4070-2004LR	4070-3004LR		
	VistaStak, 192 tips/layer, 5 layers (small size) or 3 layers (micro and large size)	4060-9024	4060-9025	4060-9026		
	Stacked Rack, 200 tips/layer, 5 layers (small size) or 3 layers (large size)		9025 ¹	9026 ¹		
	Stacked Rack, Trace Metal Certified, 200 tips/layer 5 layers (small size) or 3 layers (large size)		9022 ¹	9023 ¹		
	Box 100 tips/box (5mL) or 35 tips/box (10mL)				9048	9050
	VistaClear Box 35 tips/box					4058-6100
	VistaBulk™, 1000 tips/bag, or 250 tips/bag (5mL) or 100 tips/bag (10mL)	4058-1000	4058-2000	4058-3000	4058-5000	4058-6000
	Protectainer™ Bulk Pack, 1000 tips (small) or 750 tips (large)		4025	4026		
	Econo-Pak™ Bulk Pack, 1000 tips		4225	4226		
STERILE	VistaRak, Sterile, non-pyrogenic, certified RNase/DNase free, 192 tips/rack, 5 racks	4060-1032	4060-2132	4060-3132		
	VistaRak, Sterile, non-pyrogenic, certified RNase/DNase free, Low Retention, 192 tips/rack, 5 racks	4070-1032LR	4070-2132LR	4070-3132LR		
	VistaClear Box – Sterile, 60 tips/box (5mL) or 35 tips/box (10mL)				4058-5102	4058-6102
	VistaTip™ Individually Wrapped Sterile, 200 tips		2025	2026		
	VistaTip Individually Wrapped, Sterile, non-pyrogenic, 200 tips		2027	2028		
FILTERED	VistaRak, Filtered, Sterile, non-pyrogenic, certified RNase/DNase free, 192 tips/rack, 5 racks † 2µL tip for Ovation 0.2-2µL, Ovation 1µL and Ovation 2µL only	4060-1332 ² 4060-1333†	4060-2332	4060-3332		
	VistaRak, Filtered, Sterile, non-pyrogenic, Low-Retention, certified RNase/DNase free, 192 tips/rack, 5 racks † 2µL tip for Ovation 0.2-2µL, Ovation 1µL and Ovation 2µL only	4070-1332LR ² 4070-1333LR†	4070-2332LR	4070-3332LR		
	VistaClear Box – Filtered, Sterile, non-pyrogenic, certified RNase/DNase free, 60 tips/box (5mL) or 35 tips/box (10mL)				4058-5133	4058-6133
	VistaTip Individually Wrapped, Filtered, Sterile, non-pyrogenic, certified RNase/DNase free, 50 tips				4058-5332	4058-6332




VistaRak™



VistaStak™



Stacked Rack



VistaClear™ Box

Notes: 1. Not for use with Ovation E8 and E12 Multichannel Pipettes. 2. Not for use with Ovation QS or Ovation M .02-2µL and Ovation F2 1µL or 2µL.
3. All 5mL tips are graduated, marked at 0.1mL increments